

Recombinant Human β -Actin,ACTB (C-6His)

Catalog No: #AP75687



Package Size: #AP75687-1 10ug #AP75687-2 50ug #AP75687-3 500ug #AP75687-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------|---|
| Product Name | Recombinant Human β -Actin,ACTB (C-6His) |
| Host Species | E.coli |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | P60709 |
| Uniprot | P60709 |
| GeneID | 60; |
| Calculated MW | 42.8kDa |
| Target Sequence | <p>DDIAALVVDNGSGMCKAGFAGDDAPRAVFPSIVGRPRHQVMVGMGQKDSYVGDEAQSQRGILTLKYPYIEH GIVTNWDDMEKIWHHTFYNELRVAPEEHPVLLTEAPLNPKANREKMTQIMFETFNTPAMYVAIQAVLSLYASG RTTGIVMDSGDGVTHTVPIYEGYALPHAILRLDLAGRDLTDYLMKILTERGYSFTTTAEREIVRDIKEKLCYVALD FEQEMATAASSSSLEKSYELPDGQVITIGNERFRCPEALFQPSFLGMESCGIHETTENSIMKCDVDIRKDLNAN TVLSGGTTMYPGIADRMQKEITALAPSTMKIKIAPPKYSVWIGGSILASLSTFQQMWISKQEYDESGPSIVH RKCFLHHHHHHH</p> |
| Formulation | Supplied as a 0.2 μ m filtered solution of 10mM TrisHCl, 0.1% TritonX-100, 2mM DTT, pH 8.0. |
| Storage | Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles. |

Note: This product is for in vitro research use only