

# Recombinant Human Amyloid $\beta$ A4 Precursor Protein-Binding A3,APBA3,X11- $\gamma$ (C-6His)

Catalog No: #AP75688

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75688-1 10ug #AP75688-2 50ug #AP75688-3 500ug #AP75688-4 1mg

## Description

Product Name	Recombinant Human Amyloid $\beta$ A4 Precursor Protein-Binding A3,APBA3,X11- $\gamma$ (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O96018
Uniprot	O96018
GeneID	9546;
Calculated MW	15.48kDa
Target Sequence	MDFPTISRSPSGPPAMDLEGPRDILVPSEDLTPDSQWDPMPGGPGSLSRMELDESSLQELVQQFEALPGDLV GPSPGGAPCPLHIATGHGLASQEIADAHGLLSAEAGRDDLGLLHCEECPPSQTGPEEPLEPAPRLLLEHHHHH H
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of 20mM PB ,150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only