

Recombinant Human Cysteine Dioxygenase Type 1, CDO1 (N-6His)

Catalog No: #AP75689

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75689-1 10ug #AP75689-2 50ug #AP75689-3 500ug #AP75689-4 1mg

Description

Product Name	Recombinant Human Cysteine Dioxygenase Type 1, CDO1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q16878
Uniprot	Q16878
GeneID	1036;
Calculated MW	25.1kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMEQTEVLKPRTLADLIRILHQLFAGDEVNVEEVQAIMEAYESDPTWAMYAK FDQYRYTRNLVDQGNKFNLMILCWGEGHGSSIHDTNSHCFLKMLQGNLKETLFAWPKKSNEMVKKSER VLRENQAYINDSVGLHRVENISHTEPVSLHLYSPPFDTCHAFDQRTGHKNKVTMTFHSKFGIRTPNATSGS LENN
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mM DTT, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only