

Recombinant Human Superoxide Dismutase [Cu-Zn],SOD1,Cu-Zn SOD (N-6His)

Catalog No: #AP75690

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75690-1 10ug #AP75690-2 50ug #AP75690-3 500ug #AP75690-4 1mg

Description

Product Name	Recombinant Human Superoxide Dismutase [Cu-Zn],SOD1,Cu-Zn SOD (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00441
Uniprot	P00441
GeneID	6647;
Calculated MW	18.1kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMATKAVCVLKGDGPVQGIINFEQKESNGPVKVVWGSIKGLTEGLHGFHVHEF GDNTAGCTSAGPHFNPLSRKHGGPKDEERHVGDLGNVTADKDGADVSIEDSVISLSGDHCCIIGRTLTVHEKA DDLKGGNEESTKTGNAGSRLACGVIGIAQ
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only