

Recombinant Human MEK-Binding Protein 1,MP1,MAPKSP1 (N-6His)

Catalog No: #AP75692

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75692-1 10ug #AP75692-2 50ug #AP75692-3 500ug #AP75692-4 1mg

Description

Product Name	Recombinant Human MEK-Binding Protein 1,MP1,MAPKSP1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UHA4
Uniprot	Q9UHA4
GeneID	8649;
Calculated MW	15.78kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMADDLKRFLYKKLPSVEGLHAIWVSDRDGVPVIKVANDNAPEHALRPGFLST FALATDQGSKLGLSKNKSIIICYNTYQVVQFNRLPLVVSFIASSANTGLIVSLEKELAPLFEELRQVVEVS
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 5mM DTT, 30% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only