

Recombinant Human Growth Factor Receptor-Bound Protein 2,GRB2,ASH (C-6His)



Catalog No: #AP75704

Orders: order@signalwayantibody.com

Package Size: #AP75704-1 10ug #AP75704-2 50ug #AP75704-3 500ug #AP75704-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Growth Factor Receptor-Bound Protein 2,GRB2,ASH (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P62993
Uniprot	P62993
GeneID	2885;
Calculated MW	26.3kDa
Target Sequence	MEAIKYDFKATADDELSFKRGDILKVLNEECDQNWYKAELNGKDGFIKPNYIEMKPHPWFFGKIPRAKAEEM LSKQRHDGAFLIRESESAPGDFSLSVKFGNDVQHFKVLRDGAGKYFLWVVKFNSLNELVDYHRSTSVSRNQQ IFLRDIEQVPQQPTYVQALDFDFDPQEDGELGFRRGDFIHVMDNSDPNWWKGACHGQTGMFPRNYVTPVNRN VLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only