

# Recombinant Human Aldo-Keto Reductase 1C4,AKR1C4 (N-6His)

Catalog No: #AP75705

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75705-1 10ug #AP75705-2 50ug #AP75705-3 500ug #AP75705-4 1mg

## Description

Product Name	Recombinant Human Aldo-Keto Reductase 1C4,AKR1C4 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P17516
Uniprot	P17516
GeneID	1109;
Calculated MW	39.3kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMDPKYQRVELNDGHFMPVLGFGTYAPPEVPRNRRAVEVTKLAIEAGFRHIDS AYLYNNEEQVGLAIRSKIADGSVKREDIFYTSKLWCTFFQPQMVQPALESSLKKLQLDYVDLYLLHFPMALKPG ETPLPKDENGKVIFDTVDLSATWEVMEKCKDAGLAKSIGVSNFNRYRQLEMILNKPGLKYKPCVNQVECHPYLN QSKLLDFCKSKDIVLVAHSALGTQRHKLWDPNSPVLLEDPVLCALAKKHKRTPALIALRYQLQRGVWVLAQSY NEQRIRENIQVFEFQLTSEDMKVL DGLNRNYRYVMDFLMDHPDY PFSDEY
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only