

Recombinant Human GADD45 $\beta$ ,GADD45B (N-6His)

Catalog No: #AP75706



Package Size: #AP75706-1 10ug #AP75706-2 50ug #AP75706-3 500ug #AP75706-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GADD45 $\beta$ ,GADD45B (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O75293
Uniprot	O75293
GeneID	4616;
Calculated MW	20kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMTLEELVACDNAAQKMQTVTAAVEELLVAAQRQDRLTVGVYESAKLMNVDP DSVVLCLLAIDEEEDDIALQIHFTLIQSFCCDNDINIVRVSGMQRLAQLLGPAETQGTTEARDLHCLLVTPHT DAWKSHGLVEVASYCEESRGNNQWVPYISLQER
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only