

Recombinant Human Galactokinase, GALK1 (C-6His)

Catalog No: #AP75707



Package Size: #AP75707-1 10ug #AP75707-2 50ug #AP75707-3 500ug #AP75707-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Galactokinase, GALK1 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P51570
Uniprot	P51570
GeneID	2584;
Calculated MW	43.3kDa
Target Sequence	MAALRQPQVAELLAEARRAFREFEFGAEP ELAVSAPGRVNLIGEHTDYNQGLVLPMALELMTVLVGS PRKDGL VSLTTSEGADEPQRLQFPLPTAQRSLEPGTPRWANYVKGVIQYYP AAPLPGFSAVVVSSVPLGGGLSSSASL EVATYTFLLQQLCPDSGTIARAQVCQQA EHSFAGMPCGIMDQFISLMGQKKGHALLIDCRSLETSLVPLSDPKL AVLITNSNVRHSLASSEYPVRRRQC EEEVARALGKESLREVQLEEAARDLVSKEGFRRARHV VGEIRRTAQA AAALRRGDYRAFGRMLVESHRS LRDDYEVSCPELDQLVEAALAVPGVYGSRMTGGGFGGCTVTLLEASAAP HAMRHQEHYGGTATFYLSQAADGAKVLC LLEHHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only