

# Recombinant Human Mortality Factor 4-Like Protein 2,MORF4L2,MRGX (C-6His)



Catalog No: #AP75711

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75711-1 10ug #AP75711-2 50ug #AP75711-3 500ug #AP75711-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Mortality Factor 4-Like Protein 2,MORF4L2,MRGX (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q15014
Uniprot	Q15014
GeneID	9643;
Calculated MW	33.37kDa
Target Sequence	MSSRKQGSQPRGQQSAEEENFKKPTRSNMQRSKMRGASSGKKTAGPQQKNLEPALPGRWGGRSAENPPS GSRKTRKNKQKTPGNDDGGSTSEAPQPPRKKRARADPTVESEAFKNRMEVKVPIPEELKPWLVEDWDLV TRQQLFQLPAKKNVDAILEEYANCKKSQGNVDNKEYAVNEVVAGIKEYFNVMLGTQLLYKFERPQYAEILLAH PDAPMSQVYGAPHLRLRFVRIGAMLAYTPLDEKSLALLLGYLHDFLKYLAKNASLFTASDYKVASAEYHRKAL LEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only