

Recombinant Human ADP-Ribosylarginine Hydrolase,ADPRH,ARH1 (N-6His)



Catalog No: #AP75712

Orders: order@signalwayantibody.com

Package Size: #AP75712-1 10ug #AP75712-2 50ug #AP75712-3 500ug #AP75712-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADP-Ribosylarginine Hydrolase,ADPRH,ARH1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P54922
Uniprot	P54922
GeneID	141;
Calculated MW	41.67kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMEKYVAAMVLSAAGDALGYINGKWEFLQDGEKIHRQLAQLGGLDALDVGR WRVSDDTVMHLATAEALVEAGKAPKLTQLYYLLAKHYQDCMEDMDGRAPGGASVHNAMQLKPGKPNGWRI PFNSHEGGCGAAMRAMCIGLRFPHHSQDLTLIQVSIESGRMTHHHPTGYLGALASALFTAYAVNSRPPLQWG KGLMELLPEAKKYIVQSGYFVEENLQHWSYFQTKWENYLKLRGILDGESAPTFPESFGVKERDQFYTSLSYSG WGGSSGHDPMIAYDAVLAAGDSWKELAHRAFFHGGSDSTAAIAGCWWGVMYGFKGVSPSNYEKLEYRN RLEETARALYSLGSKEDTVISL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 50% Glycerol, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only