

# Recombinant Human Bridging Integrator 2,BIN2,BRAP1 (N-6His)

Catalog No: #AP75713

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75713-1 10ug #AP75713-2 50ug #AP75713-3 500ug #AP75713-4 1mg

## Description

Product Name	Recombinant Human Bridging Integrator 2,BIN2,BRAP1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UBW5
Uniprot	Q9UBW5
GeneID	51411;
Calculated MW	30.58kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMAEGKAGGAAGLFAKQVQKKFSRAQEKVLQKLGKAVETKDERFEQSASNF YQQQAEGHKLYKDLKNFLSAVKVMHESSKRVSETLQEIYSSEWDGHEELKAIVWNNDLLWEDYEEKLADQAV RTMEIYVAQFSEIKERIAKRGRKLVYDSARHHLEAVQNAKKKDEAKTAKAEFEFNKAQTVFEDLNQELLEELP ILYNSRIGCYVTIFQNISNLRDVFYREMSKLNHNLYEVMSKLEKQHSN
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only