

# Recombinant Human Isopentenyl Pyrophosphate Isomerase 2,,IPPI2,IDI2 (N-6His)



Catalog No: #AP75715

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75715-1 10ug #AP75715-2 50ug #AP75715-3 500ug #AP75715-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Isopentenyl Pyrophosphate Isomerase 2,,IPPI2,IDI2 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9BXS1
Uniprot	Q9BXS1
GeneID	91734;
Calculated MW	28.9kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMSDINLDWVDRRQLQRLEEMLIVVDENDKVIGADTKRNCHLNENIEKGLLHR AFSVVLFNTKNRILIQQRSDTKVTFPGYFTDSCSSHPLYNPAELEEKDAIGVRRAAQRRLQAEKGIPGEQISPED IVFMTIYHHKAKSDRIWGEHEICYLLLVRKNVTLNPDPSSETKSILYLSQEELWELLEREARGEVKVTPWLRTIAE RFLYRWWPHLDDVTPFVELHKIHRV
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mM DTT, 0.1mM PMSF, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only