

Recombinant Human Selenocysteine Lyase, SCLY (N-6His)



Catalog No: #AP75718

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75718-1 10ug #AP75718-2 50ug #AP75718-3 500ug #AP75718-4 1mg

Description

Product Name	Recombinant Human Selenocysteine Lyase, SCLY (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96I15
Uniprot	Q96I15
GeneID	51540;
Calculated MW	50.34kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMEAAVAPGRDAPAPAASQPSGCGKHNSPERKVYMDYNATTPLEPEVIQAM TKAMWEAWGNPSSPYASGRKAKDIINAARES LAKMIGGKPDIIFTSGGTESNNLVIHSVVKHFHANQTSKGH TGGHHSPVKGAKPHFITSSVEHDSIRLPLEHLVEEQVAAVTFVPVSKVSGQTEVDDILA AVRPTTRLVTIMLANN ETGIVMPVPEISQRIKALNQERVAAGLPPILVHTDAAQALGKQRVDVEDLGVDFTIVGHKIFYGPRIGALYIRGL GEFTPLYPMLFGGGQERNFRPGTENTPMIAGLGKAAELVTQNCEAYEAHMRDVRDYLEERLEAEFGQKRIHL NSQFPGTQRLPNTCNFSIRGPRLQGHVLAQCRVLMASVGAACHSDHGDQPSPVLLSYGVFPDVARNALRL SVGRSTTRAEVDLVVQDLKQAVAQLEDQA
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only