

Recombinant Human Serpin B6,Placental Thrombin Inhibitor (N-Trx, 6His)



Catalog No: #AP75719

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75719-1 10ug #AP75719-2 50ug #AP75719-3 500ug #AP75719-4 1mg

Description

Product Name	Recombinant Human Serpin B6,Placental Thrombin Inhibitor (N-Trx, 6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P35237
Uniprot	P35237
GeneID	5269;
Calculated MW	60.7kDa
Target Sequence	MSDKIIHLTDDSFDTDLKADGAILVDFWAEWCGPCKMIAPILDEIADEYQGKLTVAKLNIDQNPGTAPKYGIRGI PTLLLFKNGEVAATKVGALSKGQLKEFLDANLAGSGSGHMHSHHHSSGLVPRGSGMKETAALKFERQHMD SPDLGTDDDDKAMADIGSEFHMDVLAANGTFALNLLKTLGKDNSKNVFFSPMSMSCALAMVYMGAKGNTAA QMAQILSFNKSGGGDIHQGFQSLLEVNKTGTQYLLRVANRLFGEKSCDFLSSFRDSCQKFKYQAEEMEELDFI SAVEKSRKHINTWVAEKTEGKIAELLSPGSVDPLTRLVLVNAVYFRGNWDEQFDKENTEERLFKVSKEEKPV QMMFKQSTFKKTYIGEIFTQILVLPYVGKELNMIIMLPDETTDLRTVEKELTYEKFVEWTRLDMMDEEEVEVSLP RFKLEESYDMESVLRNLGMTDAFELGKADFSGMSQTDLSLSKVVHKSFVEVNEEGTEAAAATAAIMMMRCAR FVPRFCADHPFLFFIQHSKTNLFCGRFSSP
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only