

Recombinant Human Glutathione Synthetase, GSH Synthetase (C-6His)



Catalog No: #AP75721

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75721-1 10ug #AP75721-2 50ug #AP75721-3 500ug #AP75721-4 1mg

Description

Product Name	Recombinant Human Glutathione Synthetase, GSH Synthetase (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P48637
Uniprot	P48637
GeneID	2937;
Calculated MW	53.45kDa
Target Sequence	ATNWGSLQDKQQLEELARQAVDRALAEGLLRTSQEPTSSEVVSYAPFTLFPSLVPSALLEQAYAVQMDFNL LVDVAVSQNAAFLEQTLSSSTIKQDDFTARLFDIHKQVLKEGIAQTVFLGLNRSYMFQRSADGSPALKQIEINTIS ASFGGLASRTPAVHRHVLVLSKTKEAGKILSNPSKGLALGIAKAWELYGSPNALVLLIAQEKERNIFDQRAIE NELLARNIHVIRRTFEDISEKGSGLDQDRRLFVDGQEIADVYFRDGYMPRQYSLQNWEARLLERSHAAKCPDIA TQLAGTKKVQQELSRPGMLEMLLPQPEAVARLRATFAGLYSLDVGEEGDQAIAEALAAPSRLFVKPQREGG GNNLYGEEMVQALKQLKDSEERASYILMEKIEPEPFENCLLRPGSPARVVQCISELGIFGVYVRQEKTLMVNMK HVGHLLRTKAIEHADGGVAAGVAVLDNPYPVLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, pH 7.5 .
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only