

Recombinant Human Testin, TES (C-6His)

Catalog No: #AP75723



Package Size: #AP75723-1 10ug #AP75723-2 50ug #AP75723-3 500ug #AP75723-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Testin, TES (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UGI8
Uniprot	Q9UGI8
GeneID	26136;
Calculated MW	49.1kDa
Target Sequence	MDLENKVKKMGLGHEQGFAPCLKCKEKCEGFELHFWRKICRNCKCGQEEHDVLLSNEEDRKVKGKLFEDTK YTTLIAKLSKSDGIPMYKRNVMILTNPVAAKKNVSINTVTYEWAPPVQNLARQYMQMLPKEKQPVAGSEGAQ YRKKQLAKQLPAHDQDPSKCHELSPREVKEMEQFVKKYKSEALGVGDVKLPCEMDAQGPKQMNIPGGDRST PAAVGAMEDKSAEHKRTQYSCYCCKLSMKEGDPAIYAERAGYDKLWHPACFVCSTCHELLVDMIVFWKNEKL YCGRHYCDSEKPRCAGCDELIFSNEYTQAENQNWHLKHFFCCFDCDSILAGEIYVMVNDKPVCKPCYVKNHAV VCQGCHNAIDPEVQRVTYNNFSWHASTEFLCSCCSKCLIGQKFMPEVGMVFCVSECKKRMSLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only