

Recombinant Human Karyopherin Subunit α -2 ,KPNA2,IPOA1 (N-6His)

Catalog No: #AP75725

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75725-1 10ug #AP75725-2 50ug #AP75725-3 500ug #AP75725-4 1mg

Description

Product Name	Recombinant Human Karyopherin Subunit α -2 ,KPNA2,IPOA1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P52292
Uniprot	P52292
GeneID	3838;
Calculated MW	60.04kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMSTNENANTPAARLHRFKNKGDSTEMRRRRRIEVNVELRKAKKDDQMLKR RNVSSFPDDATSPLQENRNNQGTVNWSVDDIVKGINSSNVENQLQATQAARKLLSREKQPPIDNIIRAGLIPKF VSFLGRTDCSPIQFESAWALTNIASGTSEQTKVVVDGGAIPAFISLLASPHAHISEQAVWALGNIAGDGSVFRDL VIKYGAVDPLLALLAVPDMSSLACGYLRNLTWTLNLCRNKNPAPPIDAVEQILPTLVRLHHDDPEVLADTCW AISYLTDGPNERIGMVVKTGVVPQLVKLLGASELPVTPALRAIGNIVTGTDEQTQVVIDAGALAVFPSSLTNPKT NIQKEATWTMSNITAGRQDQIQVVNHGLVPFLVSVLSKADFKTQKEAVWAVTNYTSSGGTVEQIVYLVHCGIIE PLMNLKTAKDTKIILVILDAISNIFQAAENLGETEKL SIMIEECGGLDKIEALQNHENESVYKASLSLIEKYFSVEEE EDQNVVPETTSEGYTFQVQDQDAPGTFNF
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM TrisHCl, 1mM DTT, 20% Glycerol, pH 8.0 .
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only