

Recombinant Human D-Amino-Acid Oxidase,DAO (N-6His)

Catalog No: #AP75726

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75726-1 10ug #AP75726-2 50ug #AP75726-3 500ug #AP75726-4 1mg

Description

Product Name	Recombinant Human D-Amino-Acid Oxidase,DAO (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P14920
Uniprot	P14920
GeneID	1610;
Calculated MW	41.64kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMRVVVIGAGVIGLSTALCIHERYHSQLPLDIKVYADRFTPLTTTDVAAGLWQ PYLSDPNNPQEADWSQQTFDYLLSHVHSPNAENLGLFLISGYNLFHEAIPDPSWKDVTVLGFRKLTPRELDMFP DYGYGWFHTSLILEGKNYLQWLTERLTERGVKFFQRKVESFEEVAREGADVIVNCTGVWAGALQRDPLLQPG RGQIMKVDAPWMKHFILTHDPERGIYNSPYIIPGTQTVTLGGIFQLGNWSELNNIQDHNTIWECCRLEPTLKN ARIIGERTGFRPVRPQIRLEREQLRTGPSNTEVIHNYGHGGYGLTIHWGCALEAAKLFGRILEEKLSRMPPSH L
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM DTT, 2mM EDTA, pH 8.0 .
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only