

Recombinant Human High Mobility Group Protein B1, HMGB1



Catalog No: #AP75728

Orders: order@signalwayantibody.com

Package Size: #AP75728-1 10ug #AP75728-2 50ug #AP75728-3 500ug #AP75728-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human High Mobility Group Protein B1, HMGB1
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09429
Uniprot	P09429
GeneID	3146;
Calculated MW	10.5kDa
Target Sequence	GKGDPPKPRGKMSSYAFFVQTCREEHKKKHPDASVNFSEFSKKSERWKTMSAKEKGFEDMAKADKARY EREMKTYIPPKGETKKKF
Formulation	Lyophilized from a 0.2 µm filtered solution of 50mMHepes, 500mMNaCl, 0.5mMDTT, pH7.9 .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only