

Recombinant Human BCAS2,DAM1 (C-6His, N-T7 tag)

Catalog No: #AP75730

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75730-1 10ug #AP75730-2 50ug #AP75730-3 500ug #AP75730-4 1mg

Description

Product Name	Recombinant Human BCAS2,DAM1 (C-6His, N-T7 tag)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O75934
Uniprot	O75934
GeneID	10286;
Calculated MW	28.58kDa
Target Sequence	MASMTGGQQMGRGSMAGTGLVAGEVVVDALPYFDQGYEAPGVREAAAALVEEETRRYRPTKNYLSYLTAP DYSAFETDIMRNEFERLAARQPIELLSMKRYELPAPSSGQKNDITAWQECVNN SMAQLEHQAVRIENLELMSQ HGCNAWKVYNENLVHMIHAQKELQLRKHIQDLNWQRKNMQLTAGSKLREMESNWVSLVSKNYEIERTIVQ LENEYQIKQQHGEANKENIRQDFLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 2mM DTT, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only