

Recombinant Human Zinc Finger Protein 70,ZNF70 (N-6His)



Catalog No: #AP75732

Orders: order@signalwayantibody.com

Package Size: #AP75732-1 10ug #AP75732-2 50ug #AP75732-3 500ug #AP75732-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Zinc Finger Protein 70,ZNF70 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UC06
Uniprot	Q9UC06
GeneID	7621;
Calculated MW	53kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMEVPPATKFGETFAFENRLESQQGLFPGEDLGDPPFLQERGLEQMAVIYKEIP LGEQDEENDDYEGNFSLCSSPVQHQSIIPGTRPQDDELFGQTFLQKSDLSMCQIIHSEEPSDCDAETDRGD SGPNAPHRTPQPAKPYACRECGKAFSQSSHLLRHLVIHTGEKPYECCECGKAFSQSSHLLRHQIIHTGEKPYE CRECGKAFRQSSALTQHQKIHTGKRPYECRECGKDFSRSSSLRKHERIHTGERPYQCKECKGKSFNQSSGLS QHRKIHTLKKPHECDLGGKAFCHRSHLIRHQRIHTGKKPYKCECGKAFSQSSNLIHRKTHTGEKPYKQKQC GKAFSQSSSLIEHQRIHTGEKPYECCQCGKAFCHSSALIQHQRIHTGKKPYTCECGKAFRHRSALIEHYKTHTR EKPYVCNLCGKSFRGSSHLIRHQKIHSGEKL
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM DTT, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only