

## Recombinant Human Arginase-2, ARG2 (C-6His)

Catalog No: #AP75733



Package Size: #AP75733-1 10ug #AP75733-2 50ug #AP75733-3 500ug #AP75733-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Arginase-2, ARG2 (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P78540
Uniprot	P78540
GeneID	384;
Calculated MW	34.2kDa
Target Sequence	MHSVAVIGAPFSQGQKRKGVEHGPAAIREAGLMKRLSSLGCHLKDFGDLSTPVPKDDLYNNLIVNPRSVGLA NQELAEVVSRAVSDGYSCVTLGGDHSLAIGTISGHARHCPDLCVWVDAHADINTPLTTSSGNLHGQPVSFLL RELQDKVPQLPGFSWIKPCISSASIVYIGLRDVPPEHFILKNYDIQYFSMRDIDRLGIQKVMERTFDLLIGKRQR PIHLSFDIDAFDPTLAPATGTPVVGGLTYREGMYIAEEIHNTGLLSALDLVEVNPQLATSEEEAKTTANLAVDVIA SSFSGQTREGGHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 10mM TrisHCl, 10mM NaCl, 1mM β-mercaptoethanol, pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only