

# Recombinant Rat C-C Motif Chemokine 5,CCL5 (N-6His)

Catalog No: #AP75737

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75737-1 10ug #AP75737-2 50ug #AP75737-3 500ug #AP75737-4 1mg

## Description

Product Name	Recombinant Rat C-C Motif Chemokine 5,CCL5 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Rat
Accession No.	P50231
Uniprot	P50231
Calculated MW	10kDa
Target Sequence	MNHKVVHHHHHMDDDDKSPYGSDDTPCCFAYLSLALPRAHVKEYFYTSSKCSNLAVVVFVTRRRNRQVCANPE KKWVQEYINYLEMS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 500mM NaCl, 2mM EDTA, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only