

Recombinant Human Protein Phosphatase 1C Catalytic Subunit,PPP1CC (N, C-6His)

Catalog No: #AP75738

Orders: order@signalwayantibody.com

Package Size: #AP75738-1 10ug #AP75738-2 50ug #AP75738-3 500ug #AP75738-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Protein Phosphatase 1C Catalytic Subunit,PPP1CC (N, C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P36873
Uniprot	P36873
GeneID	5501;
Calculated MW	40.2kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMADLDKLNIDSIIQRLEVRGSKPGKNVQLQENEIRGLCLKSREIFLSQPILLEL EAPLKICGDIHGQYYDLLRLEFYGGFPPESNYLF LGDYVDRGKSLETICLLLAYKIKYPENFFLLRGNHECASI NRIYGFYDECKRRYNIKLWKTFTDCFNCLPIAAIVDEKIFCCHGGLSPDLQSMEQIRRIMRPTDVPDQGLLCDLL WSDPKDVLGWGENDRGVSFTFGAEVVAKFLHKHDLICRAHQVVEDGYEFFAKRQLVTLFSAPNYCGEFD NAGAMMSVDETLMCSFQILKPAEKKKPNATRPVTPPRGMITQAKKLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 1mM DTT, pH 8.0 .
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only