

Recombinant Human Autophagy Related 4 Homolog C,ATG4C (C-6His)



Catalog No: #AP75741

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75741-1 10ug #AP75741-2 50ug #AP75741-3 500ug #AP75741-4 1mg

Description

Product Name	Recombinant Human Autophagy Related 4 Homolog C,ATG4C (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96DT6
Uniprot	Q96DT6
GeneID	84938;
Calculated MW	53.56kDa
Target Sequence	MEATGTDEVDKLTKFISAWNNMKYSWVLKTKTYFSRNSPVLLLGKCYHFYEDDKTLPAESGCTIEDHVIA GNVEEFRKDFISRIWLTYREEFPQIEGSALTTDCGWGCTLRGTGQMLLAQGLLHFLGRAWTPDALNIENSDS ESWTSHTVKKFTASFEASLSGEREFKPTISLKETIGKYSDDHEMRNEVYHRKIISWFGDSPLALFGLHQLIEYG KKSGKKAGDWYGPVVAHILRKAVEEARHPDLQGITIYVAQDCTVYNSDVIDKQSASMTSDNADDKAVIILVPV RLGGERTNTDYLEFVKGILSLEYCVGIIGGKPKQSYFAGFQDDSLIYMDPHYCQSFVDVSIKDFPLETFHCPS PKKMSFRKMDPSCTIGFYCRNVQDFKRASEEITKMLKFSSKEKYPLFTFVNGHSRDYDFTSTTTNEEDLFSED EKKQLKRFSTEEFVLLLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 10% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only