

Recombinant Human Catalase,CAT

Catalog No: #AP75743



Package Size: #AP75743-1 10ug #AP75743-2 50ug #AP75743-3 500ug #AP75743-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Catalase,CAT
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P04040
Uniprot	P04040
GeneID	847;
Calculated MW	59.7kDa
Target Sequence	ADSRDPASDQMQRHWKEQRAAQKADVLTTGAGNPVGDKLNVITVGPRGPLLVDVVFTDEMAHFDREIPER VVHAKGAGAFGYFEVTHDITKYSKAKVFEHIGKKTPIAVRFSTVAGESGSADTVRDPGRFAVKFYTEDGNWDL VGNNTPIFFIRDPILFPSFIHSQKRNPQTHLKDPDMVWDFWSLRPESLHQVSFLFSDRGIPDGHRHMNGYGSH TFKLVNANGEAVYCKFHYKTDQGIKNSVEDAARLSQEDPDYGIRDLFNAIATGKYPSWTFYIQVMTFNQAETF PFNPFDLTKVWPHKDYPLIPVGKLVLRNPVNYFAEVEQIAFDPSNMPPGIEASPKMLQGRLFAYPDTHRHR LGPNYLHIPVNCYPYRVRVANYQRDGPMMQDNQGGAPNYYPNSFGAPEQQPSALEHSIQYSGEVRRFNTAN DDNVTQVRAFVYVNVLNEEQRKRLCENIAGHLKDAQIFIQKKAVKNFTEVHPDYGSHIQALLDKYNAEKPKNAIH TFVQSGSHLAAREKANL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only