

Recombinant Human Baculoviral IAP Repeat-Containing Protein 5,BIRC5,Survivin

Catalog No: #AP75744

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75744-1 10ug #AP75744-2 50ug #AP75744-3 500ug #AP75744-4 1mg

Description

Product Name	Recombinant Human Baculoviral IAP Repeat-Containing Protein 5,BIRC5,Survivin
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O15392
Uniprot	O15392
GeneID	332;
Calculated MW	20.4kDa
Target Sequence	MGAPTLPPAWQPFLKDHRISTFKNWPFLGCACTPERMAEAGFIHCPTENEPDLAQCFKCFKELEGWEPDDD PIEEHKKHSSGCAFLSVKKQFEELTLGEFLKLDREKAKNKIAKETNKKKKEFEETAKKVRRAIEQLAAMD
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 0.1M NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only