

Recombinant Human Microtubule-Associated Protein 1 Light Chain 3 α ,MAP1LC3A (C-6His)



Catalog No: #AP75745

Orders: order@signalwayantibody.com

Package Size: #AP75745-1 10ug #AP75745-2 50ug #AP75745-3 500ug #AP75745-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Microtubule-Associated Protein 1 Light Chain 3 α ,MAP1LC3A (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9H492
Uniprot	Q9H492
GeneID	84557;
Calculated MW	15.3kDa
Target Sequence	MPSDRPFKQRRSFADRCKEVQQIRDQHPSKIPVIIERYKGEKQLPVLDKTKFLVPDHNMSSELVKIIIRRLQLN PTQAFFLLVNQHSMSVSVSTPIADIYEQEKDEDEGFLYMVYASQETFGFLEHHHHHH
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM Tris, 20% Glycerol, 0.1M NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only