

Recombinant Human Galectin-Related Protein, LGALS1 (C-6His)

Catalog No: #AP75746

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75746-1 10ug #AP75746-2 50ug #AP75746-3 500ug #AP75746-4 1mg

Description

Product Name	Recombinant Human Galectin-Related Protein, LGALS1 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q3ZCW2
Uniprot	Q3ZCW2
GeneID	29094;
Calculated MW	20kDa
Target Sequence	AGSVADSDAVVKLDDGHLNNSLSSPVQADVYFPRILVPCGHIKGGMRPGKVKLVMGIVDLNPESFAISLTCG DSEDPPADVAVIELKAVFTDRQLLRNSCISGERGEEQSAIPYFPFIPDQPFRVEILCEHPRFRVFDGHQLDFDYH RIQTLAIDTIKINGDLQITKLGLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 100mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only