

Recombinant Human α -Parvin, PARVA, MXRA2 (C-6His)

Catalog No: #AP75747

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75747-1 10ug #AP75747-2 50ug #AP75747-3 500ug #AP75747-4 1mg

Description

Product Name	Recombinant Human α -Parvin, PARVA, MXRA2 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NVD7
Uniprot	Q9NVD7
GeneID	55742;
Calculated MW	43.2kDa
Target Sequence	MATSPQKSPSAPKSPTPKSPPSRKKDSSFLGKLGGLARRKKAKEVSELQEEGMNAINLPLSPIFELDPEDT MLEENEVRTMVDPNRSRDPKQLQELMKVLIDWINDVLVGERIIVKDLAEDLYDQQVLQKLFKLESEKLNVAEVT QSEIAQKQKLQTVLEKINETLKLPPRSIKWNVDSVHAKSLVAILHLLVALSQYFRAPIRLPDHVSIQVVVVQKREG ILQSRQIQEEITGNTEALSGRHERDAFDTLFDHAPDKLNVVKKTLITFVNKHLNKLNLEVTELETQFADGVYLVLL MGLLEGYFVPLHSFFLTPDSFEQKVLNVSFAFELMQDGGLEKPKPRPEDIVNCDLKSTLRVLYNLFTKYRNVEL EHHHHHH
Formulation	Supplied as a 0.2 μ m filtered solution of 50mM Tris, 150mM NaCl, 40% Glycerol, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only