

Recombinant Human DNA Fragmentation Factor Subunit α ,DFFA,DFF45,ICAD (N-6His)



Catalog No: #AP75748

Orders: order@signalwayantibody.com

Package Size: #AP75748-1 10ug #AP75748-2 50ug #AP75748-3 500ug #AP75748-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DNA Fragmentation Factor Subunit α ,DFFA,DFF45,ICAD (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O00273
Uniprot	O00273
GeneID	1676;
Calculated MW	38.7kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMEVTGDAGVPESGEIRTLKPCLLRRNYSREQHGVAASCLEDLRSKACDILAI DKSLTPVTLVLAEDGTIVDDDDYFLCLPSNTKFVALASNEKWAYNNSDGGTAWISQESFDVDETDGAGLKW KNVARQLKEDLSSIILLSEEDLQMLVDAPCSDLAQELRQSCATVQRLQHTLQQVLDQREEVRQSKQLLQLYLQ ALEKEGSLLSKQEESKAAFGEEVDAVDTGISRETSSDVALASHILTALREKQAPELSLSSQDLELVTKEDPKALA VALNWDIKKTETVQEACEWELALRLQQTQSLHSLRSISASKASPPGDLQNPKRARQDPT
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only