

Recombinant Human Xaa-Pro Aminopeptidase P3,XPNPEP3 (N, C-6His)

Catalog No: #AP75749

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75749-1 10ug #AP75749-2 50ug #AP75749-3 500ug #AP75749-4 1mg

Description

Product Name	Recombinant Human Xaa-Pro Aminopeptidase P3,XPNPEP3 (N, C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NQH7
Uniprot	Q9NQH7
GeneID	63929;
Calculated MW	60.2kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMPWLLSAPKLVPAVANVRGLSGCMLCSQRRYSLQVPVERRIPNRYLGQPSP FTHPHLLRPGEVTPGLSQVEYALRRHKLMSLIQKEAQQSGTDQTVVLSNPTYMSNDIPYTFHQDNNFLYL CGFQEPDSILVLQSLPGKQLPSHKAILFVPRRDPRELWDGPRSGTDGAIALTGVDEAYTLEEFQHLLPKMKA ETNMVWYDWMRPSHAQLHSDYMQPLTEAKAKSKNKVIRGVQQLIQRLRLIKSPAIEIRMQIAGKLTSAFIETM FTSKAPVEEAFLYAKFEFECRARGADILAYPPVAGGNRSNTLHYVKNNQLIKDGEMVLLDGGCESSCYVSDI TRTWPNVNGRFTAPQAELYEAVLEIQRDCLALCFPGTSLNIYSMMLTLIGQKLKDLGIMKNIKENNAFKAARKY CPHHVGHYLGMDVHDTDPMPRSLPLQPGMVITIEPGIYIPEDDKDAPEKFRGLGVRIEDDVVVTQDSPFILSAD CPKEMNDIEQICSQASLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 25mM Tris, 1mM DTT, pH 7.3.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only