

Recombinant Human Sulfotransferase 2B1,SULT2B1 (C-6His)

Catalog No: #AP75750

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75750-1 10ug #AP75750-2 50ug #AP75750-3 500ug #AP75750-4 1mg

Description

Product Name	Recombinant Human Sulfotransferase 2B1,SULT2B1 (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O00204
Uniprot	O00204
GeneID	6820;
Calculated MW	36.6kDa
Target Sequence	MDGPAEPQIPGLWDTYEDDISEISQKLPGEYFRYKGVPPVGLYSLESISLAENTQDVRDDDIFIITYPKSGTTW MIEIICLILKEGDPSWIRSVPIWERAPWCETIVGAFSLPDQYSPRLMSSHLPIQIFTKAFFSSKAKVIYMGRNPRD VVVSLYHYSKIAGQLKDPGTPDQFLRDFLKGVEVQFGSWFDHIKGWLRMKGKDNFLFITYEELQQDLQGSVERI CGFLGRPLGKEALGSVVAHSTFSAMKANTMSNYTLLPPSLLDHRRGAFLRKGVCGDWKNHFTVAQSEAFDR AYRKQMRGMPTFPWDEV DHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 0.1M NaCl, 1mM EDTA, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only