

Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase D,PPID,PPlase D (N, C-6His)



Catalog No: #AP75752

Orders: order@signalwayantibody.com

Package Size: #AP75752-1 10ug #AP75752-2 50ug #AP75752-3 500ug #AP75752-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase D,PPID,PPlase D (N, C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q08752
Uniprot	Q08752
GeneID	5481;
Calculated MW	43.9kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMSHPSPQAKPSNPSNPRVFFDVIDIGGERVGRIVLELFADIVPKTAENFRALC TGEKGIGHTTGKPLHFKGCPFHRIIKKFMIQGGDFSNQNGTGGESIYGEKFEDENFHYKHDREGLLSMANAGR NTNGSQFFITTVPTPHLDGKHVVFQVIKIGIVARILENVEVKGEKPAKLCVIAECGELKEGDDGGIFPKDGSG DSHPDFPEDADIDLKDVKILLITEDLKNIGNTFFKSQNWEMAIIKKYAEVLRVYDSSKAVIETADRKLQPIALSC VLNIGACKLKMSNWQGAIDSCLEALELDPSTNKALYRRAQGWQGLKEYDQALADLKKAAQGIAPEDKAIQAELL KVKQKIKAQKDKKAVYAKMFAVEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only