

Recombinant Human α -Crystallin A Chain,CRYAA (C-6His)

Catalog No: #AP75755

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75755-1 10ug #AP75755-2 50ug #AP75755-3 500ug #AP75755-4 1mg

Description

Product Name	Recombinant Human α -Crystallin A Chain,CRYAA (C-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P02489
Uniprot	P02489
GeneID	102724652;1409;
Calculated MW	20.9kDa
Target Sequence	MDVTIQHPWFKRTLGPFPYSRLFDQFFGEGLEFYDLLPFLSSTISPYRQSLFRTVLDSGISEVRSRDRDKFVIFL DVKHFSPEDLTVKVQDDFVEIHGKHNERQDDHGYISREFHRRYRLPSNVDQSALSCSLSDGMLTFCGPKIQT GLDATHAERAIPVSREEKPTSAPSSLEHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, 2mM EDTA, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only