

# Recombinant Human Methionine Aminopeptidase 1D, MetAP1D, MAP1D (N, C-6His)



Catalog No: #AP75757

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75757-1 10ug #AP75757-2 50ug #AP75757-3 500ug #AP75757-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Methionine Aminopeptidase 1D, MetAP1D, MAP1D (N, C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q6UB28
Uniprot	Q6UB28
GeneID	254042;
Calculated MW	35.4kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMRQRDISHSIVLPAAVSSAHPVPKHIKKPDYVTTGIVPDWGDSIEVKNEDQIQ GLHQACQLARHVLLLAGKSLKVDMTTEEIDALVHREIISHNAYPSPLGYGGFPKSVCTSVNNVLCHGIPDSRPL QDGDIIINIDVTVYYNGYHGDTSETFLVGNVDECCKLVEVARRCRDEAIAACRAGAPFSVIGNTISHITHQNGF QVCPHFVGHGIGSYFHGHPEIWHHANDSDLPMEEGMAFTIEPIITEGSPEFKVLEDAWTVVSLDNQRSAQFEH TVLITSRGAQILTKLPHEALEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 50mM Tris, 100mM NaCl, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only