

# Recombinant Human FGF1 Intracellular-Binding Protein, FIBP



Catalog No: #AP75759

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75759-1 10ug #AP75759-2 50ug #AP75759-3 500ug #AP75759-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human FGF1 Intracellular-Binding Protein, FIBP
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O43427
Uniprot	O43427
GeneID	9158;
Calculated MW	41.8kDa
Target Sequence	AMADIGSMTSELDIFVGNTTLLIDEDVYRLWLDGYSVTDVAVALRVRSGILEQTGATAAVLQSDTMDHYRTFHMLE RLLHAPPKLLHQLIFQIPPSRQALLIERYAFDEAFVREVLGKKLSKGTCKDLDDISTKTGITLTKSCRRQFDNFKR VFKVVEEMRGLVDNIQQHFLLSDRLARDYAAIVFFANNRFETGKKKLQYLSFGDFAFCAELMIQNWTLGAVD SQMDDMDMDLDKEFLQDLKELKVLVADKDLLLHKSLVCTALRGKLGVFSEMEANFNLSRGLVNVAACKLTH NKDVRDLFVDLVEKFVEPCRSDHWPLSDVRFFLNQYSASVHSLDGFRRHQALWDRYMGTLRGCLLRLYHD
Formulation	Supplied as a 0.2 µm filtered solution of 50mM TrisHClpH8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only