

# Recombinant Human Deoxycytidine Kinase,DCK (N-6His, T7 tag)

Catalog No: #AP75760

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75760-1 10ug #AP75760-2 50ug #AP75760-3 500ug #AP75760-4 1mg

## Description

Product Name	Recombinant Human Deoxycytidine Kinase,DCK (N-6His, T7 tag)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P27707
Uniprot	P27707
GeneID	1633;
Calculated MW	34kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMASMTGGQQMGRGSMATPPKRSCPSFSASSEGTRIKKISIEGNIAAGKSTF VNILKQLCEDWEVPEPVARWCNVQSTQDEFEELTMSQKNGGNLQMMYEKPERWSFTFQTYACLSRIRAQ LASLNGKDKDAEKPVLFFERSVYSDRYIFASNLYESECMNETEWTIYQDWHDFWMNNQFGQSLELDGIYLQAT PETCLHRIYLRGRNEEQGIPLEYLEKLHYKHESWLLHRTLKTNFDYLQEVPIILTDVNEFDKDKYESLVEKVKEF LSTL
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHClpH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only