

Recombinant Human D- β -Hydroxybutyrate Dehydrogenase Mitochondrial, BDH1 (N-6His)



Catalog No: #AP75762

Orders: order@signalwayantibody.com

Package Size: #AP75762-1 10ug #AP75762-2 50ug #AP75762-3 500ug #AP75762-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human D- β -Hydroxybutyrate Dehydrogenase Mitochondrial, BDH1 (N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q02338
Uniprot	Q02338
GeneID	622;
Calculated MW	40.3kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMLATRLSRPLSRLPGKTLSDRENGARRPLLLGSTSFIPIGRRTYASAAEPV GSKAVLVTGCDSGFGFALAKHLHSGFLVFAGCLMKDKGHDGVKELDSLNSDRLRTVQLNVCSSSEVEKVV IVRSSLKDPEKGMWGLVNNAGISTFGEVEFTSLETYKQVAEVLNLTGTVRMTKSFLPLIRRAKGRVVNISSMLG RMANPARSPYCITKFGVEAFSDCLRYEMYPLGVKVSVEPGNFIAATSLYSPESIQAIKKMWEELPEVVRKDY GKKYFDEKIAKMETYCSSGSTDTSVIDAVTHALTATTPYTRYHPMDYYWLRMQIMTHLPGAISDMIYIR
Formulation	Supplied as a 0.2 μ m filtered solution of 50mM Tris, 1mM DTT, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only