

Recombinant Human Follistatin-Related Protein 1, FSTL1 (C-6His, E. coli)

Catalog No: #AP75775

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75775-1 10ug #AP75775-2 50ug #AP75775-3 500ug #AP75775-4 1mg

Description

Product Name	Recombinant Human Follistatin-Related Protein 1, FSTL1 (C-6His, E. coli)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q12841
Uniprot	Q12841
GeneID	11167;
Calculated MW	33.8kDa
Target Sequence	MEEELRSKSKICANVFCGAGRECAVTEKGEPTCLCIEQCKPHKRPVCGSNGKTYLNHCELHRDACTGSKIQV DYDGHCKEKKSVSPSASPVCYQSNRDELRRRIIQWLEAEIIPDGWFSKGSNYSEILDKYFKNFDNGDSRLDS SEFLKFVEQNETAINITTPDQENKLLRGLCVDALIELSDENADWKLSFQEFKCLNPSFNPPEKKCALEDETY ADGAETEVDNRCVCACGNWVCTAMTCDGKNQKGAQTQTEEMTRYVQELQKHQETAECTKRVSTKEIVD HHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only