

Recombinant Human Ubiquitin-Conjugating Enzyme E2 I,UBE2I (N-GST)



Catalog No: #AP75780

Orders: order@signalwayantibody.com

Package Size: #AP75780-1 10ug #AP75780-2 50ug #AP75780-3 500ug #AP75780-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ubiquitin-Conjugating Enzyme E2 I,UBE2I (N-GST)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P63279
Uniprot	P63279
GeneID	7329;
Calculated MW	44.4kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKLTQSMAIIRYI ADKHNMLGGCPKERAEISMLEGAVLDIRYGVSRAYSDFETLKVDFLSKLPPEMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIQIDKYLKSSKYIAWPLQGWQATFGGGDHPKSDL VPRGSHMSGIALSRLAQERKAWRKDHPFGFVAVPTKNPDGTMNLMNWECAIPGKKGTPWEGGLFKLRMLFK DDYPSSPPKCKFEPPLFHPNVYPSGTVCLSILEEDKDWRPATIKQILLGIQELLNEPNIQDPAQAEAYTIYCQNR VEYEKRVRAQAKKFAPS
Formulation	Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only