

Recombinant Human dUTP Pyrophosphatase,dUTPase

Catalog No: #AP75782

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75782-1 10ug #AP75782-2 50ug #AP75782-3 500ug #AP75782-4 1mg

Description

Product Name	Recombinant Human dUTP Pyrophosphatase,dUTPase
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P33316-2
Uniprot	P33316
GeneID	1854;
Calculated MW	17.7kDa
Target Sequence	MPCSEETPAISPSKRARPAEVGGMQLRFARLSEHATAPTRGSARAAGYDLYSAYDYTIPMEKAVVKTDIQIAL PSGCYGRVAPRSGLAAKHFIDVGAGVIDEDYRGNVGVVLFNFGKEKFEVKKGDRIAQLICERIFYPEIEEVQAL DDTERGSGGFGSTGKN
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only