

TPH1 (Phospho-Ser58) Polyclonal Antibody

Catalog No: #12328

Package Size: #12328-1 50ul #12328-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TPH1 (Phospho-Ser58) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-TPH1 (S58) Polyclonal Antibody detects endogenous levels of TPH1 protein only when phosphorylated at S58.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human TPH1 around the phosphorylation site of S58.
Target Name	TPH1
Modification	Phospho
Other Names	TPH1; TPH; TPRH; TRPH; Tryptophan 5-hydroxylase 1; Tryptophan 5-monoxygenase 1
Accession No.	Swiss-Prot: P17752NCBI Gene ID: 7166
Uniprot	P17752
GeneID	7166;
Target Species	human
SDS-PAGE MW	51kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300

ELISA: 1/10000

Not yet tested in other applications

Images



Western Blot analysis of 293 cells using Phospho-TPH1 (S58) Polyclonal Antibody

Note: This product is for in vitro research use only