

Recombinant Human Ubiquitin-Conjugating Enzyme E2 B,UBE2B,HR6B



Catalog No: #AP75804

Orders: order@signalwayantibody.com

Package Size: #AP75804-1 10ug #AP75804-2 50ug #AP75804-3 500ug #AP75804-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ubiquitin-Conjugating Enzyme E2 B,UBE2B,HR6B
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P63146
Uniprot	P63146
GeneID	7320;
Calculated MW	17.5kDa
Target Sequence	GSHMSTPARRRLMRDFKRLQEDPPVGVSGAPSENNIMQWNAVIFGPEGTPFEDGTFKLVIEFSEEYPNKPPT VRFLSKMFHPNVYADGSICLDILQNRWSPTYDVSSILTSIQSLLDEPNPNPANSQAAQLYQENKREYEKRVSA IVEQSWNDS
Formulation	Supplied as a 0.2 µm filtered solution of 10mM HEPES, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only