

Recombinant Human Cystathionine γ -Lyase, CTH

Catalog No: #AP75805



Package Size: #AP75805-1 10ug #AP75805-2 50ug #AP75805-3 500ug #AP75805-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Cystathionine γ -Lyase, CTH
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P32929
Uniprot	P32929
GeneID	1491;
Calculated MW	44.7kDa
Target Sequence	GHMQEKDASSQGFLPHFQHFATQAIHVGQDPEQWTSRAVVPISLSTTFKQGAPGQHSGFEYSRSGNPTRN CLEKAVAALDGAKEYCLAFASGLAATVTITHLLKAGDQIICMDDVYGGTNRVFRQVASEFGLKISFVDCSKIKLLE AAITPETKLVWIETPTNPTQKVIDIEGCAHIVHKHGDILVVDNTFMSPYFQRPLALGADISMYSATKYMNGHSDV VMGLVSVNCESLHNRLRFLQNSLGAVSPIDCYLCNRGLKTLHVRMEKHFKNGMAVAQFLESNPWVEKVIYP GLPSHPQHELVKRQCTGCTGMVTFYIKGTLQHAEIFLKNLKLFTLAESLGGFESLAELPAIMTHASVLKNDRDV LGISDTLIRLSVGLLEDEEDLLEDLDQALKAHPPSGSHS
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, 20% Glycerol, pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only