

Recombinant Human Insulin-Like Growth Factor II,IGF2



Catalog No: #AP75806

Orders: order@signalwayantibody.com

Package Size: #AP75806-1 10ug #AP75806-2 50ug #AP75806-3 500ug #AP75806-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Insulin-Like Growth Factor II,IGF2
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01344
Uniprot	P01344
GeneID	3481;
Calculated MW	8.91kDa
Target Sequence	MFPAMPLSSLFVNAYRPSETLCGGELVDTLQFVCGDRGFYFSRPASRVSRRSRGIVECCFRSCDLALLETY CATPAKSE
Formulation	Lyophilized from a 0.2 µm filtered solution of 5mM Hac, pH ~3.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only