

# Recombinant Human Death Domain-Containing Protein CRADD,CRADD

Catalog No: #AP75809

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75809-1 10ug #AP75809-2 50ug #AP75809-3 500ug #AP75809-4 1mg

## Description

Product Name	Recombinant Human Death Domain-Containing Protein CRADD,CRADD
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P78560
Uniprot	P78560
GeneID	8738;
Calculated MW	23kDa
Target Sequence	GSHMEARDKQVLRSLRLELGAEVLVEGLVLQYLYQEGILTENHIQEINAQTTGLRKTMLLLDILPSRGPKAFDTF LDSLQEFPPWREKLKKAREEAMTDLPAGDRLTGIPSHILNSSPSDRQINQLAQRLGPEWEPMVLSLGLSQTDI YRCKANHPHNVSQVVEAFIRWRQRFQKQATFQSLHNGLRAVEVDPSLLLHMLE
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only