

Recombinant Human Astrocytic Phosphoprotein PEA-15,PEA15

Catalog No: #AP75810

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75810-1 10ug #AP75810-2 50ug #AP75810-3 500ug #AP75810-4 1mg

Description

Product Name	Recombinant Human Astrocytic Phosphoprotein PEA-15,PEA15
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q15121
Uniprot	Q15121
GeneID	8682;
Calculated MW	15.3kDa
Target Sequence	GSHMAEYGTLLQDLTNNITLEDLEQLKSACKEDIPSEKSEEITGSAWFSFLESHNKLDKDNLSYIEHIFEISRRP DLLTMVVDYRTRVLKISEEDELDTKLTRIPSAKKYKDIIRQPSEEEIILAPPPKKA
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only