

Recombinant Human  $\gamma$ -Synuclein, SNCG

Catalog No: #AP75812



Package Size: #AP75812-1 10ug #AP75812-2 50ug #AP75812-3 500ug #AP75812-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human $\gamma$ -Synuclein, SNCG
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O76070
Uniprot	O76070
GeneID	6623;
Calculated MW	13.6kDa
Target Sequence	GSHMDVFKKGFSLAKEGVVGAVEKTKQGVTEAAEKTKEGVMYVGAKTKENVVQSVTSVAEKTKEQANAVSEA VVSSVNTVATKTVEEAENIAVTSGVVRKEDLRPSAPQQEASKEKEEVAEEAQSGGD
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only